



MedCare
SPA PREMIUM PRODUCTS

Innovation
for your
Health!

Watercleaning
Skin Care
Intim hygienie
Pulmonology

CONTENTS

WATERCLEANING

SPA KDF/Ag Carbon Filter	CODE 3372
SPA Ag Carbon Filter	CODE 3389
SPA Innovative KDF Filter	CODE 3396
SPA Innovative Ag Carbon Filter	CODE 3303
SPA FILTER CAPSULE	CODE 3365
SPA FAUCET	CODE 3860
SPA Wall Water Filter Faucet	CODE 3884
SPA Spout for Faucet Water Pipe	CODE 3303
SPA Massage Shower head	CODE 3426
SPA Spray Shower Head	CODE 3433
SPA Multifunction Shower Head Neige	CODE 3181
SPA Multifunction Shower Head Minerve	CODE 3174
SPA Multifunction Shower Head Clear	CODE 3136

SKIN CARE

SPA Derma capsule	CODE 3341
SPA Multifunction Shower Head Neige	CODE 3181
SPA Multifunction Shower Head Minerve	CODE 3174
SPA Multifunction Shower Head Clear	CODE 3136

INTIM HYGIENE

SPA Douche	CODE 3907
SPA DOUCHE Applicator	CODE 3235
SPA Complex Salt capsula	CODE 3969
SPA Lactic Acid Capsule	CODE 3914
SPA Calendula and Salvia Capsule	CODE 3945
SPA Magistral Capsule	CODE 3952
SPA Multifunction Shower Head Neige	CODE 3181
SPA Multifunction Shower Head Minerve	CODE 3174
SPA Multifunction Shower Head Clear	CODE 3136

PULMONOLOGY

Pulmotrainer	CODE 3105
--------------	-----------

SPA KDF / AG CARBON FILTER

The capsule is made of medical PE and the cartridge contains 50% KDF 55[®] PROCESS MEDIUM and 50% silver coated activated carbon granules. These two special materials ensure perfect filtration in applications where you need to filter and purify water contaminated by heavy metals and various types of bacteria. It filters and kills bacteria and fungi. When the KDF cartridge of the capsule turns from the original golden colour into dark brown or black, the capsule must be replaced.

The filter capsules may only be used for potable tap water!

SPA AG CARBON FILTER

The capsule is made of medical PE and the cartridge contains 100% silver coated activated carbon granules. It is recommended for water filtration where unwanted odors, flavors and other organic matters need to be removed from tap water. It filters and kills bacteria and fungi. When the colour of the granules turns from greyish black into yellowish grey, the capsule must be replaced. **The filter capsules may only be used for drinkable tap water!**



CODE 3372



CODE 3389

WATERCLEANING



SPA INNOVATIVE KDF FILTER

The capsule is made of 100% PA-lined brass, which provides such a large filtration surface that it multiplies the efficiency and life expectancy of conventional KDF 55[®] PROCESS MEDIUM filters. The cartridge of the capsule contains 100% KDF 55[®] PROCESS MEDIUM. In this case filtration aims to remove chlorine and heavy metals. When the colour of the brass turns brownish black, the capsule must be replaced. **The filter capsules may only be used for potable tap water!**

CODE 3396



SPA INNOVATIVE AG CARBON FILTER

The capsule is made of 100% brass, the cartridge contains 100% silver coated activate carbon granules. In addition to being an efficient filter (eliminating unwanted flavors and odors), it also kills bacteria and fungi, preserving the useful minerals and taste of water. When the color of the brass turns brownish black, the capsule must be replaced. **The filter capsules may only be used for drinkable tap water!**

WATERCLEANING

SPA FILTER CAPSULE

The capsule filters cold and hot drinking water in both the bathroom and the kitchen. Can be used for up to 4 months depending on water quality and frequency of use. It is recommended that you check the capsule periodically – when its contents go from their original gold color to dark brown or black the capsule needs to be replaced. During use the SPA Filter Capsule turns the chlorine found in water into chloride, filters out heavy metals and other contaminants (lead, chrome, mercury, arsenic, iron etc.) and kills fungi, bacteria, therefore protecting your skin from dryness and irritation. Filled with a special KDF 55[®] PROCESS MEDIUM filter material.



CODE 3365

SPA FAUCET

The SPA FAUCET combined with the SPA FILTER CAPSULE is a mini water cleaning device in hospitals, in laboratories, a industrial setting and at home.

Material: copper with chrome plating



CODE 3860

WATERCLEANING



SPA WALL WATER FILTER FAUCET

The SPA Faucet is made up of two pieces: the base of the faucet and the spout piece of the device. The Spa Filter capsule can be inserted into the device by removing the spout piece.

Material: copper with chrome plating

CODE 3884



SPA SPOUT FOR FAUCET WATER PIPE

Water saving aerator and spout head that can be turned 360°, which makes everyday use easier. There is no more need for buying bottled water - with SPA you save and protect the environment.

The spout can be purchased separately, it can be connected to any other wall mixer with a ½-inch thread.

Material: copper with chrome plating

CODE 3853

WATERCLEANING

SPA MASSAGE SHOWER HEAD

The Spa Massage Water Cleaning Shower is an elegant premium category shower head made of solid copper, which is known to have antibacterial qualities.

Included in the device is a replaceable pre-filter that provides easy access for cleaning and lengthens the lifetime of the capsules used with the shower head.

Material: copper with chrome plating

SPA SPRAY SHOWER HEAD

The Spa Spray Shower Head also includes a valve found only in premium product lines that prevents water from streaming backwards. For instance, in case of reduced water pressure, the built in valve stops water from flowing backwards, therefore preventing potential bacterial inflow and spreading.

Material: copper with chrome plating



CODE 3426



CODE 3433

WATERCLEANING



SPA SHOWER HEADS WITH SWAROWSKI CRISTALS

The Spa Multifunction Shower comes in an elegant design and is primarily used for water filtering. The specially designed, hollow handle and the detachable shower head make up the body of the product. Similarly to most shower heads the Spa Multifunction Shower Head has to be attached to the shower hose. The device is compatible with most standard shower hoses and **ONLY** works with the Spa capsules.

Only available for special order.

WATERCLEANING

SPA MULTIFUNCTION SHOWER HEADS

The SPA Multifunction Shower Heads are suitable for one of the following functions:

- Watercleaning
- Skin care
- Intimate hygiene

Except every day using the Spa Multifunction Shower Heads can also provide treatment options with the use of **dermatological, vaginal, anal** Spa capsules.



Neige

CODE 3181



Minerve

CODE 3174



Clear

CODE 3136

SKIN CARE



MULTIFUNCTION SHOWER HEADS WITH SPA DERMA CAPSULE

The Spa Derma Capsule is specifically designed skin moisturizing and dermatological treatment capsules that can be placed into the hollow handle of the SPA Multifunction Shower Head.

The Spa Derma capsule can be filled with physician prescribed medicinal treatment. The capsule gradually and continuously dispenses its contents during shower, ensuring even distribution on the surface of your skin and providing a simple yet effective treatment.

SPA DERMA CAPSULE

The Spa Derma capsule is available for custom dermatologic treatments. Medications prescribed by your physician can be placed into the empty capsule. The capsule is designed to gradually dissolve its contents into the water and dispense the treatment drop by drop. The medicated rinse covers can be applied to your entire body enabling you to complete the treatment quickly, efficiently, on your own.

The Spa Derma capsule can **ONLY** be used with the Spa Multifunction Shower head and the Spa Shower. The capsule can **ONLY** be filled with physician prescribed medicinal treatment.



CODE 3341

INTIM HYGIENE



SPA DOUCHE (SPA Multifunction Shower + SPA Applicator)

SPA Multifunction Shower with Spa Douche Applicator and the supplementing Spa capsules make everyday hygiene along with gynecological and proctological treatments possible.

Can be attached to the bath or shower faucet and is compatible with most models worldwide.

The treatments can be done in the peaceful surroundings of your own home, without any need for external help.

The Spa Douche ONLY works with Spa capsules (sold separately).

CODE 3907

SPA DOUCHE APPLICATOR

The Spa Douche Applicator is an optional accessory (sold separately) to the Spa Multifunction Shower and serves as replacement applicator head for the Spa Douche.


The Spa Douche Applicator does not work with any other devices.

The Spa Douche Applicator turns the Spa Multifunction Shower into a unique tool used for gynecological and anal treatments in the comfort of your home. The applicator is an ergonomically designed instrument, made with tissue-friendly, soft plastic.



CODE 3235

INTIM HYGIENE



ENEMA

Lack of regular exercise and poor diet leads to constipation in 70% of our population. Sitting for long periods of time daily along with other causes can result in a lack of bowel movement for days or painful passing of the stool, which can lead to anal tears and harm to the epidermis of the anus. These injuries heal slowly, can get infected, and can be accompanied by hemorrhoids. These painful symptoms can be avoided with the occasional administrations of enemas that introduce liquid through the anus, resulting in rectal stimulation and bowel movement.

The Spa Douche provides a safe, simple solution that enables you to administer enemas in the comfort of your home. Other devices required the help of a person with rectal irrigation experience, which has prevented a lot of people from using enemas. However, now our one of a kind technology lets you complete this process without requiring assistance.

ANAL TREATMENT WITH ONE-PHASE CAPSULE

SPA COMPLEX CAPSULE

The Spa Complex Capsule is specifically designed to be used with the SPA MULTIFUNCTION SHOWER WITH APPLICATOR or the SPA DOUCHE.

The NaCl content of the capsule dissolves in water, letting 8-10 dl liquid into your anus, after which the device stops water flow automatically. The risk of high water pressure entering your anus is eliminated thanks to the water draining channels built into the device. Maintain good anal hygiene during everyday cleansing.

•Recommended uses:

- Enema, Anal irrigation
- Constipation
- Improving hemorrhoid symptoms
- Before and after anal intercourse



CODE 3969



VAGINAL TREATMENT WITH TWO-PHASE CAPSULES

The SPA Capsule's top part is filled with NaCl salt crystals. These salt crystals dissolve with the inflow of approximately 0.5 dl water during the first phase of the treatment and clean the walls of the vagina or anus. After the first phase is complete we insert the second phase contents located in the bottom part of the capsule into the now clean vagina or anus. This entire procedure lasts no more than mere seconds and the device shuts off automatically upon completion of the treatment.

For vaginal use:

- Lactic Acid Capsule
- Calendula and Salvia Capsule
- Magistral

LACTIC ACID CAPSULE

FOR VAGINAL USE

The SPA Lactic Acid Capsule is a natural treatment for restoring and maintaining a healthy vaginal flora and pH balance. Acidic environment prevents bacterial infections, fungus, and the growth of other parasites and organisms in the vagina.

Recommended uses:

- Vaginal discharge
- Infections caused by sexual intercourse
- Fungal infections caused by swimming pools
- Abrasions
- Vaginal discharge caused by antibiotic treatments
- Vaginal discomfort caused by hormonal changes
- Restoring vaginal pH balance
- Preventing infections



CODE 3914

INTIM HYGIENE



CODE 3945

CALENDULA AND SALVIA CAPSULE

FOR VAGINAL USE

Clinical studies indicate that calendula and salvia successfully treat many unpleasant smelling vaginal discharge types that are not caused by microorganisms such as fungi and bacteria. Calendula and salvia are believed to improve vaginal resistance, reduce inflammation, and help heal micro damage of the vagina as well. The lactic acid content of these herbs also aids in restoring and maintaining a healthy vaginal pH balance.

Recommended uses:

- Vaginal discharge
- Infections caused by sexual intercourse
- Fungal infections caused by swimming pools
- Abrasions
- Vaginal discharge caused by antibiotic treatments
- Vaginal discomfort caused by hormonal changes
- Restoring vaginal pH balance

MAGISTRAL CAPSULE

FOR VAGINAL OR ANAL USE

The SPA Magistral Capsule is for localized vaginal and anal treatment. Can ONLY be used with physician prescribed ingredients! The Spa Magistral capsule can ONLY be used with liquid, gel, or soft cream contents prescribed by your physician.❏

Recommended uses:

- Custom treatment
- Safe, exact dosage
- Liquid, gel, cream

Two-phase capsule that gives you the option to customize your treatment. The capsule's top part is filled with NaCl salt crystals and the bottom part can be filled with physician prescribed ingredients. Please consult your doctor for proper usage information.



CODE 3952

WHY IS BETTER THE SPA DOUCHE THAN OTHER PRODUCTS?

The Spa Douche has a built in pressure regulator:

The vaginal wall lies against the head of the applicator. The existing products in the market have a round, plain surface that does not have a drain to let the water out of the vagina. These are gravity fed instruments and approximately 0.7-1 Bar pressure by hand (or 3.5-4 Bar pressure with the other device, which can be connected to the water system) forces the water into the vagina and it cannot exit due to the muscle surrounding the instrument. As a result, the vagina first expands like a balloon then the water finds its way up towards the uterus. This is very dangerous because it takes the infected vaginal fluids that we are trying to flush out even deeper into the body. It can often lead to pelvic inflammatory disease and other infections as well. The Spa Douche has deep ridges that drain the water out of the vagina to ensure that the water leaves immediately. In case of anal usage the same ridges prevent overpressure and guarantee the exit of excess water.

The Spa Douche guarantees one-way water flow

:If the head of the hand operated, balloon instruments not removed from the vagina after the pressure is applied a suction effect will come into play. Upon release the infected fluid enters the balloon immediately and the already infected fluids are pushed back into the body. With the Intimate Shower the water flows in one direction through the capsules, which allows only the filtered and cleaned water to enter the vagina. A suction effect will not happen even when the water pressure seizes

The Spa Douche only rinses in a sideways direction:

Other instruments let the water enter the vagina through a hole located on the head of the applicator and shoot the rinsing liquid directly into the vagina with great pressure. On the head of the Spa Douche applicator the holes involved in the rinsing process only spray water sideways. The rinsing fluid only reaches the sides of the vagina, then leaves through the above mentioned ridges immediately.

The Spa Douche regulates water intake:

The pH of the vagina should not be washed out completely, which is why it is important that we don't use more than 2 ounces of water for the rinsing process. For both the manually operated instrument and the attached device the user can decide how much water to irrigate. The moment the contents of the capsule dissolve is not obvious, therefore, it is possible that the user could be rinsing out the herbs and minerals with excess water. The Two-phase Spa capsules are engineered to dissolve with 2 ounces of water. Once the capsules are empty the Spa Douche is designed to stop automatically without any intervention from the user. Therefore, over rinsing is not possible. (Capsules used for enema are designed to dissolve with 4 ounces of water. Once the capsule is empty water flow stops automatically.)

The Spa Douche provides a controlled amount of medication and consistent care quality:

If the patient is responsible for preparing and portioning out the ingredients for the treatment it is possible that not the proper amount of medication is dispersed with the device. This can not happen with the Spa Douche because every capsule is filled with a controlled amount of ingredients and the quality of the treatment is consistent. The only exception is the Magistral capsule, in which only the top part is pre-filled by the manufacturer. The bottom part should be filled by physician prescribed ingredients in quantities allowed by the capsule's space limitations. The capsules are disposable and cannot be refilled.

The Spa Douche ensures that drug intake is in the upper part of the vagina:

The most effective way to apply medication to the vagina is through localized treatment with the use of liquids. However, this is difficult to do. Using brushes to apply medication in the vagina can only be done by physicians. At home only cream and suppository can be used which is unnecessary and contains harmful ingredients. Applying these is complicated and can result in medicating the wrong place. The Spa Douche ensures drug intake into the upper part of the vagina in the liquid form. Active agents from the bottom part of the capsule flow through the inner channel of the applicator head and into the vagina freely.

PULMONOLOGY



Pulmonary function improver device for asthmatic and respiratory problems Postoperative use lung training to increase vital capacity and lung function, avoiding postoperative complication (for example postoperative inflammation).

Compared to previous models, this device is not sensitive to motions during usage. There is no need to adjust the setting angle for proper use. The needed amount of exhaled air and the right pressure can be adjusted during use. As a result of exhalation a vibration causes oscillations inside the house. This causes a vibration in the airway walls, which in turn facilitates the movement of secretion up the airways, consequently improving lung function and oxygenation.

PULMOTRAINER

Recommended for patients with the following diagnosis:

- Chronic Obstructiv Pulmonary Disease (COPD)
- Chronic bronchitis
- Mucoviscidosis
(inherited disorder that involves mucus and sweat glands)
- Bronchiectasia
(Abnormal permanent dilation of the bronchial tubes)
- Emphysema
- Asthma bronchiale

For the following symptoms:

- Chronic cough
- Productive cough

Material: plastic



CODE 3105

Distributor:

Nemes Medical AB.

352 34 Vaxjö, SWEDEN

Tradgardsgatan 16C

gynec@hotmail.com

info@medcarespa.com

www.medcarespa.com

www.pulmotrainer.se

T: + 46 727 410092

Manufacturer:

MedCare Ltd.

1116 Budapest,

Építész str. 8-12.

HUNGARY

T: + 36 30 650 3520

info@medcareltd.com

www.medcareltd.com